

DATA SHEET

Selekt[™] Off Feed

Presentation A beige powder, each 785 grammes of which has the following

composition.

| Calcium propionate | 500 g |
|------------------------------|-------|
| Magnesium sulphate 7-hydrate | 9 |
| Potassium chloride | |
| Dry yeast | _ |
| Gentian root | • |

Total calcium content: 108 g.

Uses For inappetant cows. Also provides rehydration.

Instructions for feeding

Selekt[™] Off Feed should be dissolved in warm water at a rate of one 785-gramme sachet in 20 litres. The solution should be delivered into the rumen with a Selekt[™] pump-drencher. If additional rehydration is required, give a solution of Selekt[™] Restore in 20 litres of warm water as well, i.e. 40 litres in total.

Contra-indications and warnings

Do not use less than 20 litres of water to dissolve the contents of the sachet. Do not give more Selekt™ Off Feed until 12 hours have elapsed.

nave claps

Further information Selekt™ Off Feed is to be used in accordance with veterinary

advice.

Legal category Complementary feeding stuff.

Package quantity Box containing 12 x 785 grammes

Marketed by Nimrod Veterinary Products Ltd

25a, High Street Moreton-in-Marsh

GL56 0AF

United Kingdom

Telephone: +44 1608 652593

e-mail: nimrod@nimrodvet.co.uk website: www.nimrodvet.co.uk

V6-2016-09